



The Controlled Medical Devices ( CMD ) Inspection is to be completed at the **beginning of each shift.**

### CMD Inspection Report

Check for the **presence and condition** of all components.

For the Month of:	Year:
Apparatus No:	
Defibrillator Serial No:	
EpiPen 0.3 mg ( Yellow Cap ) Lot No:	
EpiPen 0.15mg ( Green Cap ) Lot No:	

Acceptable

Damaged / Missing / Unacceptable

Zoll AEDPro Defibrillator																																
Platoon on Duty																																
Date	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
CONTENTS	Case																															
	Defibrillator																															
	Cables ( Intact )																															
	Unit Battery																															
	Spare Battery																															
	Electrode Pad Kits (2)																															
	→ Sealed																															
	→ Not expired																															
Self Test Indicator OK																																
Initials																																
Epinephrine Auto Injector																																
CONTENTS	EPI Pen (Yellow) 0.3 mg																															
	EPI Pen (Green) 0.15 mg																															
	Not Expired																															
	Medication Clear *																															
	No Condensation Evident																															
	No Discolouration/ Chunks Evident																															

**Note:**

- \* The viewing window is naturally translucent on the EpiPen® & EpiPen Jr.®
- Upon completion of daily CMD inspection, the officer is required to document “CMD Checks Complete” within the comment section of Quatro
- Store active record in an area readily accessible by all Fire Station personnel
- Store completed record in the fire station for one year after completion or until requested by TFS Mechanical Division, EMS Section personnel or Sunnybrook Centre for Prehospital Medicine